



## **ANOGEN** - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

### **Human IgG Datasheet**

<b>Catalogue No</b>	IG30001
<b>Product name</b>	Human Immunoglobulin G (IgG)
<b>Description</b>	Human Immunoglobulin G (IgG) is prepared from normal human serum and purified by Protein A affinity column
<b>Formulation</b>	Lyophilized from a solution of 0.01M PBS, pH7.4
<b>Reconstitution</b>	Double distilled water is recommended to adjust the final concentration to 1.00mg/ml
<b>Applications</b>	Research and further manufacturing
<b>Storage</b>	Store at -20° C
<b>References</b>	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

**This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administration to humans and animals for use in diagnostic and therapeutic procedures.**

**Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited**