



## **ANOGEN** - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

### **Mouse IgG Datasheet**

|                       |   |
|-----------------------|---|
| <b>Catalogue No</b>   | IG41002   |
| <b>Product name</b>   | Mouse Immunoglobulin G (IgG)  |
| <b>Description</b>    | Mouse Immunoglobulin G (IgG) is prepared from normal mouse serum and purified by Protein A affinity column  |
| <b>Formulation</b>    | Lyophilized from a solution of 0.01M PBS, pH7.4   |
| <b>Reconstitution</b> | Double distilled water is recommended to adjust the final concentration to 1.00mg/ml  |
| <b>Applications</b>   | Research and further manufacturing  |
| <b>Storage</b>        | Store at -20° C   |
| <b>References</b>     | If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude. |

**This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administration to humans and animals for use in diagnostic and therapeutic procedures.**

**Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited**