



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Mouse IgG2b Datasheet

Catalogue No	IG41004
Product name	Mouse Immunoglobulin G2b (IgG2b)
Description	Mouse Immunoglobulin G2b (IgG2b) is prepared from normal mouse serum and purified by Protein A affinity column
Formulation	Lyophilized from a solution of 0.01M PBS, pH7.4
Reconstitution	Double distilled water is recommended to adjust the final concentration to 1.00mg/ml
Applications	Research and further manufacturing
Storage	Store at -20° C
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administration to humans and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited