



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Rat IgG Datasheet

Catalogue No	IG42003
Product name	Rat Immunoglobulin G (IgG)
Description	Rat Immunoglobulin G (IgG) is prepared from normal rat serum and purified by Protein A affinity column
Formulation	Lyophilized from a solution of 0.01M PBS, pH7.4
Reconstitution	Double distilled water is recommended to adjust the final concentration to 1.00mg/ml
Applications	Research and further manufacturing
Storage	Store at -20° C
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administration to humans and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited