



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

HRP-conjugated Human Triiodothyronine (T3) Datasheet

Catalogue No	MO-M40044T
Product name	HRP conjugated Human Triiodothyronine (T3)
Product Type	Conjugated Antigen
Description	Human Triiodothyronine (T3), HRP conjugated.
Target Molecule	Antibody of Triiodothyronine (T3)
Specificity	Reactive to Antibody of Triiodothyronine (T3)
Reactivity	Human, others not tested.
Source	Triiodo-L-thyronine sodium salt [Sigma Cat No: T2752]
Formulation	Liquid form with 50% PBS and 50% glycerol.
Applications	ELISA, this product can be used as tracer in competence ELISA with our anti-Triiodothyronine monoclonal antibody, [Clone T3-1; Cat No: MO-M40044B]. The suggested tracer concentration is 1 : 1000 dilution or follow the document with shipping,
Research Area	Metabolism, Hormone.
Storage	Store at -20° C. Aliquot and store at -20°C; avoid thaw and freeze cyclically.
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited

Rev: S7.5(01)