



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Human Chorionic Gonadotrophin (hCG) Datasheet

Catalogue No	AP30002
Product name	Human chorionic gonadotrophin (hCG)
Description	HCG from human urine
Source	Human urine
Format	Lyophilized powder
Formulation	Lyophilized powder
Reconstitution	Double distilled water is recommended to adjust the final concentration to 1.00 mg/mL
Stability	Should be stored below 0° C. Reconstituted form is best stored refrigerated at 4° C
Activity	Biopotency: 3200IU/mg. Immunopotency: 6000IU/mg
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administration to humans and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited