



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

CRP Antigen Datasheet

Catalogue No	AP30030
Product name	CRP antigen
Description	Highly purified CRP antigen from human fluid
Source	Human
Format	Liquid.
Concentration	2.75 mg/mL
Formulation	Liquid form: colorless, clear, no offensive odour. 0.1M Tris 0.2M NaCl pH 7.5 +/- 0.2 0.1% NaN ₃ 2mm CaCl ₂
Storage Conditions	Optimum Short Term Storage Temperature: 4° C Optimum Long Term Storage Temperature: -20° C
Applications	This material is offered for research, laboratory or further manufacturing purposes Hazard Testing: HIV1, HIV2 Results: All donor units tested negative for HIV 1,2 HbsAg Results: All donor units tested negative for HbsAg HCV Results: All donor units tested negative for HCV
Precautions	No test guarantees a product to be non-infectious. Therefore, all material derived from human fluids or tissues should be considered as potentially infectious
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administration to humans and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited