



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Mouse anti-Human TNF- α Monoclonal Antibody [C4C7] Datasheet

Catalogue No	MO-C40003H
Product name	Mouse anti-Human TNF- α Monoclonal Antibody [C4C7]
Clone No	C4C7
Product type	Primary antibodies
Description	Mouse monoclonal antibody to human Tumor Necrosis Factor alpha (TNF- α), Protein G affinity purified
Target protein	Human TNF- α
Immunogen	Purified recombinant human TNF- α
Specificity	This antibody recognizes recombinant human TNF- α in ELISA and Western Blotting.
Reactivity	Human, other species not tested
Clonality	Monoclonal
Source	Mouse
Myeloma	Sp2/0-Ag14
Subclass	IgG1
Light Chain	Kappa
Formulation	Lyophilized from a solution in 0.01M PBS, pH 7.0
Reconstitution	Double distilled water is recommended to adjust the final concentration to 1.00 mg/ml
Applications	ELISA and Western Blotting.
Research	Growth Factors, Hormones and their Receptors.
Storage	Store at -20° C
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.


Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited

S7.5(01)



ANOGEN - A Division of YES Biotech Laboratories Ltd.
2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Figure: Western blot analysis of Human TNF-alpha using anti-human TNF-alpha monoclonal antibody



This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited

S7.5(01)