



## **ANOGEN - A Division of YES Biotech Laboratories Ltd.**

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

### **Mouse anti-Human EGF Monoclonal Antibody [S-146] Datasheet**

<b>Catalogue No</b>	MO-C40044F
<b>Product name</b>	Mouse anti-human EGF monoclonal antibody [S-146]
<b>Clone No</b>	S-146
<b>Product type</b>	Primary antibodies
<b>Description</b>	Mouse monoclonal antibody to human epidermal growth factor (EGF), Protein G affinity purified
<b>Target protein</b>	Human EGF
<b>Immunogen</b>	Purified recombinant human EGF, 6kDa, 53 amino acid residues
<b>Specificity</b>	Reactive with natural and recombinant human EGF
<b>Reactivity</b>	Human, others not tested
<b>Cross-reactivity</b>	Does not cross react with other human cytokines tested such as IL-1 $\beta$ , IL-8, MCAF, TGF- $\beta$ and SAA
<b>Clonality</b>	Monoclonal
<b>Source</b>	Mouse
<b>Myeloma</b>	Sp2/0-Ag14
<b>Subclass</b>	IgG1
<b>Light Chain</b>	Kappa
<b>Formulation</b>	Lyophilized from a solution in 0.01M PBS, pH 7.2
<b>Reconstitution</b>	Double distilled water is recommended to adjust the final concentration to 1.00 mg/ml
<b>Applications</b>	ELISA
<b>Matched Pair</b>	Clone No.: S-134    Cat. No.: MO-C40044D Clone No.: S-146    Cat. No.: MO-C40044F
<b>Research</b>	Growth Factor / Cytokine
<b>Storage</b>	Store at -20° C
<b>References</b>	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

**This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.**

**Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited**