




**ANOGEN** - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

### Mouse anti-Human bFGF Monoclonal Antibody [F-74] Datasheet

<b>Catalogue No</b>	MO-C40046B
<b>Product name</b>	Mouse anti-human bFGF monoclonal antibody [F-74]
<b>Clone No</b>	F-74
<b>Product type</b>	Primary antibodies
<b>Description</b>	Mouse monoclonal antibody to human fibroblast growth factor basic (bFGF), protein G affinity purified
<b>Target protein</b>	Human bFGF
<b>Immunogen</b>	Purified recombinant human bFGF
<b>Reactivity</b>	Human, others not tested
<b>Clonality</b>	Monoclonal
<b>Source</b>	Mouse
<b>Subclass</b>	IgG2a
<b>Formulation</b>	Lyophilized from a solution in 0.01M PBS, pH 7.2
<b>Reconstitution</b>	Double distilled water is recommended to adjust the final concentration to 1.00 mg/ml
<b>Applications</b>	ELISA, Western Blot
<b>Research</b>	Growth Factor
<b>Storage</b>	Store at -20° C
<b>References</b>	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

**Figure:** Western blot analysis of Human bFGF using anti-human bFGF monoclonal antibody



**This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.**

**Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited**

S7.5(01)