




ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Mouse anti-Human bFGF Monoclonal Antibody [F-198] Datasheet

Catalogue No	MO-C400461
Product name	Mouse anti-human bFGF monoclonal antibody [F-198]
Clone No	F-198
Product type	Primary antibodies
Description	Mouse monoclonal antibody to human fibroblast growth factor basic (bFGF), protein G affinity purified
Target protein	Human bFGF
Immunogen	Purified recombinant human bFGF
Reactivity	Human, others not tested
Clonality	Monoclonal
Source	Mouse
Subclass	IgG1
Formulation	Lyophilized from a solution in 0.01M PBS, pH 7.2
Reconstitution	Double distilled water is recommended to adjust the final concentration to 1.00 mg/ml
Applications	ELISA, Western Blot
Research	Growth Factor
Storage	Store at -20° C
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

Figure: Western blot analysis of Human bFGF using anti-human bFGF monoclonal antibody



This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited

S7.5(01)