



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Mouse anti-Human M-CSF Monoclonal Antibody [116] Datasheet

Catalogue No	MO-C40048A
Product name	Mouse anti-human M-CSF monoclonal antibody [116]
Clone No	116
Product type	Primary antibodies
Description	Mouse monoclonal antibody to human macrophage colony stimulating factor (M-CSF), Protein G affinity purified
Target protein	Human M-CSF
Immunogen	Purified recombinant human M-CSF (rhM-CSF)
Specificity	This monoclonal antibody is specifically reactive to rhM-CSF
Reactivity	Human, others not tested
Cross-reactivity	Does not cross react with BSA or other human cytokines tested such as IL-1 β , IL-8, IL-16, EGF, G-CSFR, GM-CSF, MCP-1, MCP-3, TGF- β and TNF- α
Clonality	Monoclonal
Source	Mouse
Myeloma	Sp2/0-Ag14
Subclass	IgG2b
Light Chain	Kappa
Formulation	Lyophilized from a solution in 0.01M PBS, pH 7.2
Reconstitution	Double distilled water is recommended to adjust the final concentration to 1.00 mg/ml
Applications	ELISA
Research	Growth Factor
Storage	Store at -20 $^{\circ}$ C
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administration to humans and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited