



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Mouse anti-HIV-1 P24 (Core Protein) Monoclonal Antibody [1027] Datasheet

Catalogue No	MO-I40002J
Product name	Mouse anti-HIV-1 P24 (core protein) monoclonal antibody [1027]
Clone No	1027
Product type	Primary antibodies
Description	Mouse monoclonal antibody (MAb) to HIV-1 P24 (core protein), Protein G affinity purified
Target protein	Human HIV-1 P24 (core protein)
Immunogen	Purified recombinant HIV-1 P24 (core protein)
Specificity	This MAb reacts with HIV-1 P24 (core protein)
Reactivity	Human, others not tested
Cross-reactivity	No cross activity with HIV-1 GP41 and HIV-2 GP36
Clonality	Monoclonal
Source	Mouse
Myeloma	Sp2/0-Ag14
Subclass	IgG1
Light Chain	Kappa
Formulation	Lyophilized from a solution in 0.01M PBS, pH 7.2
Reconstitution	Double distilled water is recommended to adjust the final concentration to 1.00 mg/ml
Applications	ELISA Western Blot: 0.1-0.5µg/mL of antibody used with 0.1µg/lane of antigen in reduced SDS-PAGE gel then blotted on PVDF membrane.
Research	Virology
Storage	Store at -20° C
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited

S7.5(01)