



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Mouse anti-HEV Monoclonal Antibody [S-13] Datasheet

Catalogue No	MO-I40019B
Product name	Mouse anti-HEV monoclonal antibody [S-13]
Clone No	S-13
Product type	Primary antibodies
Description	Mouse monoclonal antibody (MAb) to human Hepatitis E Virus, Protein G affinity purified
Target protein	Human hepatitis E virus (HEV) antigen
Immunogen	A highly antigenic polypeptide consisting of essential sequences at least 30 residues in length which were selected from gene encoding ORF3 region of a Chinese HEV strain
Reactivity	Human, others not tested
Cross-reactivity	No reactivity with HCV-S, HCV-NS3, HCV-NS4
Clonality	Monoclonal
Source	Mouse
Myeloma	Sp2/0-Ag14
Subclass	IgM
Light Chain	Kappa
Formulation	Lyophilized
Concentration	Double distilled water is recommended to adjust the final concentration to 1.00 mg/ml
Applications	ELISA
Research	Virology
Storage	Store at -20° C
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administration to humans and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited