



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Mouse anti-Human CD11 α Monoclonal Antibody [A11] Datasheet

Catalogue No	MO-L40028
Product name	Mouse anti-human CD11 α monoclonal antibody [A11]
Clone No	A11
Product type	Primary antibodies.
Description	Mouse monoclonal antibody to human CD11 α , Protein G affinity purified.
Target protein	Human CD11 α
Immunogen	Human cytolytic T lymphocytes.
Specificity	This antibody is reactive to human CD11 α .
Reactivity	Human, others not tested.
Clonality	Monoclonal
Source	Mouse
Subclass	IgG1
Light Chain	Kappa
Formulation	Lyophilized from a solution in 0.01M PBS, pH 7.0
Reconstitution	Double distilled water is recommended and to adjust the final concentration to 1.00 mg/ml
Applications	EIA, Immunohistochemistry.
Research	Immunology, cell type.
Storage	Store at -20° C
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited

S7.5(01)