



## **ANOGEN - A Division of YES Biotech Laboratories Ltd.**

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

### **Mouse anti-Human Transferrin Receptor Monoclonal Antibody [CY-TFR] Datasheet**

<b>Catalogue No</b>	MO-40042
<b>Product name</b>	Mouse Monoclonal Antibody to Human Transferrin Receptor[CY-TFR]
<b>Clone No</b>	CY-TFR
<b>Product type</b>	Primary antibodies
<b>Description</b>	Mouse monoclonal antibody to human Transferrin Receptor, Protein G affinity purified
<b>Target protein</b>	Human Transferrin Receptor
<b>Immunogen</b>	Purified human Transferrin Receptor
<b>Specificity</b>	Reactive with human Transferrin Receptor
<b>Reactivity</b>	Human, others not tested
<b>Cross-reactivity</b>	Not available
<b>Clonality</b>	Monoclonal
<b>Source</b>	Mouse
<b>Subclass</b>	IgG1
<b>Formulation</b>	Lyophilized from a solution in 0.01M PBS, pH 7.2
<b>Reconstitution</b>	Double distilled water is recommended and to adjust the final concentration to 1.00 mg/ml
<b>Applications</b>	ELISA.
<b>Research</b>	Oncology and other related area.
<b>Storage</b>	Store at -20° C
<b>References</b>	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

**This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.**

**Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited**