



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Mouse anti-Human SAA Monoclonal Antibody [607] Datasheet

Catalogue No	MO-C40028E
Product name	Mouse anti-human SAA monoclonal antibody [607]
Clone No	607
Product type	Primary antibodies
Description	Mouse monoclonal antibody to human serum amyloid A (SAA), Protein G affinity purified
Target protein	Human SAA
Immunogen	Highly Purified recombinant human SAA (MW: 12kDa)
Specificity	Reacts with natural and recombinant human SAA
Reactivity	Human, others not tested
Cross-reactivity	Does not show any cross-reaction with other human cytokines or growth factors tested such as IL-1 β , IL-8, MCAF, TGF- β and EGF
Clonality	Monoclonal
Source	Mouse
Myeloma	Sp2/0-Ag14
Subclass	IgG2 α
Light Chain	Kappa
Formulation	Lyophilized from a solution in 0.01M PBS, pH 7.2
Reconstitution	Double distilled water is recommended and to adjust the final concentration to 1.00 mg/ml
Applications	ELISA, Western Blot: An antibody concentration of 0.1-0.5 μ g will allow visualization of 100 ng/lane of recombinant human SAA
Matched Pair	Clone 607, Cat No.: MO-C40028E Clone 115, Cat No.: MO-C40028A
Research	Biochemistry
Storage	Store at -20 $^{\circ}$ C
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited