



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Mouse anti-hCG-a Monoclonal Antibody [CH7] Datasheet

Catalogue No	MO-M40010B
Product name	Mouse anti-hCG-a monoclonal antibody [CH7]
Clone No	CH7
Product type	Primary antibodies
Description	Mouse monoclonal antibody to human chorionic gonadotropin-alpha (hCG-a), Protein G affinity purified
Target protein	Human chorionic gonadotropin-alpha (hCG-a)
Immunogen	Purified human chorionic gonadotropin
Specificity	Reactive to chorionic gonadotropin-alpha (hCG-a)
Reactivity	Human, others not tested
Clonality	Monoclonal
Source	Mouse
Subclass	IgG1
Formulation	Lyophilized from a solution in 0.01M PBS, pH 7.2
Reconstitution	Double distilled water is recommended and to adjust the final concentration to 1.00 mg/mL
Applications	ELISA, Radio-immunoassay, Purification of HCG by Affinity Chromatography, Qualitative and Quantitative Immuno-precipitation
Research	Hormone
Storage	Store at -20° C
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administration to humans and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited