



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Mouse anti-Human beta 2 Micro Globulin Monoclonal Antibody (D6E4) Data Sheet

Catalogue No	MO-M40021B
Product name	Mouse anti-Human beta 2 Micro Globulin (MG) Monoclonal Antibody (D6E4)
Clone No	D6E4
Product type	Primary antibodies
Description	Mouse anti-Human beta 2 (MG) Monoclonal Antibody), Protein G affinity purified
Target protein	Human Human beta 2 (MG)
Immunogen	Recombinant Human beta 2 (MG)
Reactivity	Human, others not tested
Clonality	Monoclonal
Source	Mouse
Subclass	IgG1
Light Chain	Kappa
Formulation	Lyophilized from a solution in 0.01M PBS, pH 7.2
Reconstitution	Double distilled water is recommended and to adjust the final concentration to 1.00 mg/mL
Research Applications	Biochemistry ELISA
Storage	Store at -20° C
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited

S7.5(01)