



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Mouse anti-Human CRP Monoclonal Antibody [S5G1] Datasheet

Catalogue No	MO-M40023E
Product name	Mouse anti-human CRP monoclonal antibody [S5G1]
Clone No	S5G1
Product type	Monoclonal antibody
Description	Mouse monoclonal antibody to human C Reactive Protein (CRP), Protein G affinity purified
Target protein	Human CRP
Immunogen	Purified CRP from human pleural fluid
Specificity	Reactive with human CRP in immunoblot and ELISA
Reactivity	Human, others not tested
Clonality	Monoclonal
Source	Mouse
Subclass	IgG1
Formulation	Lyophilized from a solution in 0.01M PBS, pH 7.2
Reconstitution	Double distilled water is recommended and to adjust the final concentration to 1.00 mg/mL
Applications	ELISA
Research	Biochemistry, Inflammation
Storage	Store at -20° C
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administration to humans and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited