



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Mouse anti-Human T3 Monoclonal Antibody [S-191] Datasheet

Catalogue No	MO-M40030A
Product name	Mouse anti-human Triiodothyronine (T3) monoclonal antibody [S-191]
Clone No	S-191
Product type	Primary antibodies
Description	Mouse monoclonal antibody to human T3, Protein G affinity purified
Target protein	Human triiodothyronine (T3)
Immunogen	T3- BSA conjugate
Reactivity	Human, others not tested
Clonality	Monoclonal
Source	Mouse
Myeloma	Sp2/0-Ag14
Subclass	IgG2a
Light Chain	Kappa
Formulation	Lyophilized from a solution in 0.01M PBS, pH 7.2
Reconstitution	Double distilled water is recommended and to adjust the final concentration to 1.00 mg/mL
Applications	ELISA
Research	Hormone
Storage	Store at -20° C
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administration to humans and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited