



**ANOGEN** - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

**Mouse anti-Human Factor IX Monoclonal Antibody [F9-1] Datasheet**

<b>Catalogue No</b>	MO-M40037
<b>Product name</b>	Mouse anti-human Factor IX monoclonal antibody [F9-1]
<b>Clone No</b>	F9-1
<b>Product type</b>	Primary antibodies
<b>Description</b>	Mouse monoclonal antibody to human Factor IX, Protein G affinity purified
<b>Target protein</b>	Human Factor IX
<b>Immunogen</b>	Human Factor IX
<b>Reactivity</b>	Human, others not tested
<b>Clonality</b>	Monoclonal
<b>Source</b>	Mouse
<b>Subclass</b>	IgG1
<b>Light Chain</b>	Kappa
<b>Formulation</b>	Lyophilized from a solution in 0.01M PBS, pH 7.2
<b>Concentration</b>	Double distilled water is recommended and to adjust the final concentration to 1.00 mg/mL
<b>Application</b>	ELISA
<b>Storage</b>	Store at -20° C
<b>References</b>	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

**This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administration to humans and animals for use in diagnostic and therapeutic procedures.**

**Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited**