



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Mouse anti-Human PTP1-C Monoclonal Antibody [PTY11] Datasheet

Catalogue No	MO-M40043B
Product name	Mouse anti-human PTP1-C monoclonal antibody [PTY11]
Clone No	PTY11
Product type	Primary antibodies
Description	Mouse monoclonal antibody to human protein-tyrosine phosphatase 1 C (PTP1-C), Protein G affinity purified
Target protein	Human PTP1-C
Immunogen	Purified recombinant protein corresponding to full-length PTP1-C
Specificity	This antibody recognizes single 67 kDa PTP1-C in human cells that express the enzyme
Reactivity	Human, others not tested
Clonality	Monoclonal
Source	Mouse
Myeloma	Sp2/0-Ag14
Subclass	IgG1
Light Chain	Kappa
Formulation	Lyophilized from a solution in 0.01M PBS, pH 7.2
Reconstitution	Double distilled water is recommended and to adjust the final concentration to 1.00 mg/mL
Application	ELISA, Western Blot, Immunocytochemistry
Research	Biochemistry
Storage	Store at -20° C
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administration to humans and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited