



**ANOGEN** - A Division of YES Biotech Laboratories Ltd.  
2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

### Mouse anti-Human CD11 $\alpha$ Monoclonal Antibody [A11] Datasheet

<b>Catalogue No</b>	MO-M40051
<b>Product name</b>	Mouse anti-human CD11 $\alpha$ monoclonal antibody [A11]
<b>Clone No</b>	A11
<b>Product type</b>	Primary antibodies.
<b>Description</b>	Mouse monoclonal antibody to human CD11 $\alpha$ , Protein G affinity purified.
<b>Target protein</b>	Human CD11 $\alpha$
<b>Immunogen</b>	Human cytolytic T lymphocytes.
<b>Specificity</b>	This antibody is reactive to human CD11 $\alpha$ .
<b>Reactivity</b>	Human, others not tested.
<b>Clonality</b>	Monoclonal
<b>Source</b>	Mouse
<b>Subclass</b>	IgG1
<b>Light Chain</b>	Kappa
<b>Formulation</b>	Lyophilized from a solution in 0.01M PBS, pH 7.0
<b>Reconstitution</b>	Double distilled water is recommended and to adjust the final concentration to 1.00 mg/ml
<b>Applications</b>	EIA, Immunohistochemistry.
<b>Research</b>	Immunology, cell type.
<b>Storage</b>	Store at -20° C
<b>References</b>	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

**This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.**

**Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited**

S7.5(01)