



## **ANOGEN - A Division of YES Biotech Laboratories Ltd.**

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

### **Mouse anti-Human NSE(145) Monoclonal Antibody Datasheet**

<b>Catalogue No</b>	MO-M40092A
<b>Product name</b>	Mouse anti-human Neuron Specific Enolase (NSE) monoclonal antibody.
<b>Clones</b>	145
<b>Product type</b>	Primary antibody.
<b>Description</b>	Mouse monoclonal antibody to NSE, protein G affinity purified.
<b>Target protein</b>	Human NSE
<b>Immunogen</b>	Recombinant human NSE
<b>Specificity</b>	This monoclonal antibody reacts with recombinant human NSE
<b>Reactivity</b>	Human, others not tested
<b>Clonality</b>	Monoclonal
<b>Source</b>	Mouse
<b>Subclass</b>	IgG2b
<b>Light chain</b>	Kappa
<b>Formulation</b>	Lyophilized from a solution in 0.01M PBS, pH 7.2
<b>Applications</b>	ELISA and Western Blot
<b>Reconstitution</b>	Double distilled water is recommended and to adjust the final concentration to 1 mg/ml.
<b>Storage</b>	Store at -20° C. After reconstitution, aliquot and store in -20° C, avoid thaw and freeze cycles.
<b>Research</b>	Nero-Study
<b>References</b>	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

**This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.**

**Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited**

Rev:0607