



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Mouse anti-Human NSE(1G6) Monoclonal Antibody Datasheet

Catalogue No	MO-M40092B
Product name	Mouse anti-human Neuron Specific Enolase (NSE) monoclonal antibody.
Clones	1G6
Product type	Primary antibody.
Description	Mouse monoclonal antibody to NSE, protein G affinity purified.
Target protein	Human NSE
Immunogen	Recombinant human NSE
Specificity	This monoclonal antibody reacts with recombinant human NSE
Reactivity	Human, others not tested
Clonality	Monoclonal
Source	Mouse
Subclass	IgM
Light chain	Kappa
Formulation	Lyophilized from a solution in 0.01M PBS, pH 7.2
Applications	ELISA and Western Blot
Storage	Store at -20° C, after reconstitution, aliquot and store in -20° C, avoid thaw and freeze cycles.
Reconstitution	Double distilled water is recommended and to adjust the final concentration to 1 mg/ml
Research	Nero-study
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited

Rev:0607