



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Mouse anti- human beta amyloid peptide N-terminal antibody [NT 5B8] Datasheet

Catalogue No	MO-M40093C
Product name	Mouse anti-human beta amyloid peptide N-terminal antibody[NT 5B8]
Clone No	NT 5B8
Product type	Primary antibodies
Description	Mouse monoclonal antibody (MAb) against the N-terminal sequence of human beta amyloid peptides, Protein G affinity purified
Target protein	human beta amyloid peptides
Immunogen	Human beta amyloid peptide 1-42
Specificity	This antibody recognizes the N-terminal sequence (DAEFRHDS) of human beta amyloid peptide 40 and peptide 42.
Reactivity	Human
Clonality	Monoclonal
Source	Mouse
Isotype	IgG1, k
Formulation	Lyophilized from a solution in 0.01M PBS pH7.2
Reconstitution	Double distilled water is recommended and to adjust the final concentration to 1.00 mg/ml
Applications	ELISA, Immuno-blotting
Research	Aging and neurodegenerative diseases.
Storage	Store at -20° C
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.


Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited

S7.5(01)



ANOGEN - A Division of YES Biotech Laboratories Ltd.
2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Western blot-beta amyloid 1-42 N terminal antibody (5B8)



1: Abeta40 peptide , 2: Abeta42 peptide

ascites (5B8) at 1:5000 dilution

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited

S7.5(01)