



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Mouse anti- human beta amyloid peptide N-terminal antibody [NT 1H3] Datasheet

| | |
|-----------------------|---|
| Catalogue No | MO-M40094B |
| Product name | Mouse anti-human beta amyloid peptide N-terminal antibody[NT 1H3] |
| Clone No | NT 1H3 |
| Product type | Primary antibodies |
| Description | Mouse monoclonal antibody (MAb) against the N-terminal sequence of human beta amyloid peptides, Protein G affinity purified |
| Target protein | human beta amyloid peptides |
| Immunogen | Full length human beta amyloid peptide 43. |
| Specificity | This antibody recognizes the N-terminal sequence of human beta amyloid peptides. |
| Reactivity | Human |
| Clonality | Monoclonal |
| Source | Mouse |
| Isotype | IgG1, k |
| Formulation | Lyophilized from a solution in 0.01M PBS pH7.2 |
| Reconstitution | Double distilled water is recommended and to adjust the final concentration to 1.00 mg/ml. Avoid repeated freeze and thaw cycle. |
| Applications | Immunohistochemistry |
| Research | Aging and neurodegenerative diseases. |
| Storage | Store at -20° C |
| References | If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude. |

Figure 1 : IHC experiment using NT-1H3 on the brain tissue from transgenic Alzheimer's disease mouse model. The plaques showed strong positive signal and the background was clear.


This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited

S7.5(01)



ANOGEN - A Division of YES Biotech Laboratories Ltd.
2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023



NT 1H3
(1:100)

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited

S7.5(01)