



## **ANOGEN - A Division of YES Biotech Laboratories Ltd.**

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

### **Mouse anti- human beta amyloid peptide 40 antibody [CV9 7B10] Datasheet**

<b>Catalogue No</b>	MO-M40095A
<b>Product name</b>	Mouse anti-human beta amyloid peptide 40 antibody[CV9 7B10]
<b>Clone No</b>	CV9 7B10
<b>Product type</b>	Primary antibody, detection antibody of matched pair for ELISA
<b>Description</b>	Mouse monoclonal antibody (MAb) against the C-terminal sequence of human beta amyloid peptide 40, Protein G affinity purified
<b>Target protein</b>	Human beta amyloid peptide 40
<b>Immunogen</b>	The C-terminal amino acid sequence (GLMVGGVV) of human beta amyloid peptide 40, conjugated with KLH.
<b>Specificity</b>	This antibody recognizes the C-terminal peptide (GLMVGGVV) of human beta amyloid peptide 40, and full length beta amyloid peptide 40.
<b>Reactivity</b>	Human
<b>Clonality</b>	Monoclonal
<b>Source</b>	Mouse
<b>Isotype</b>	IgG1, k
<b>Formulation</b>	Lyophilized from a solution in 0.01M PBS pH7.2
<b>Reconstitution</b>	Double distilled water is recommended to adjust the final concentration to 1.00 mg/ml. Avoid repeated freeze and thaw cycle.
<b>Applications</b>	ELISA
<b>Research</b>	Aging and neurodegenerative diseases.
<b>Storage</b>	Store at -20° C
<b>References</b>	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

**This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.**

**Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited**

S7.5(01)