



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Mouse anti-Human Neuron Specific Enolase (NSE) Monoclonal Antibody Datasheet[171]

Catalogue No	MO-M40091C
Product name	Mouse anti-human Neuron Specific Enolase monoclonal antibody 171
Clone No	171
Product type	Primary antibodies.
Description	Mouse monoclonal antibody (MAb) to human Neuron Specific Enolase (NSE), Protein G affinity purified
Target protein	Human Neuron Specific Enolase (NSE)
Immunogen	Recombinant Human Neuron Specific Enolase (NSE)
Specificity	Reacts with recombinant human Neuron Specific Enolase (NSE)
Reactivity	Human, others not tested
Cross-reactivity	Not available
Clonality	Monoclonal
Source	Mouse
Subclass	IgG1, k
Formulation	Lyophilized from a solution in 0.01M PBS, pH 7.2
Reconstitution	Double distilled water is recommended and to adjust the final concentration to 1.00 mg/ml
Applications	ELISA
Research	Neurology
Storage	Store at -20° C
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited

S7.5(01)