



ANOGEN - A Division of YES Biotech Laboratories Ltd.
2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023


Mouse anti-HA Tag Monoclonal Antibody [HA.C5] Datasheet

Catalogue No	MO-H40004A
Product name	Mouse anti-HA Tag Monoclonal Antibody [HA.C5]
Clone No	HA.C5
Product type	Purified primary antibodies
Description	Mouse monoclonal antibody against HA tag, an epitope of influenza hemagglutinin.
Target protein	HA tagged recombinant protein
Immunogen	A synthetic peptide (YPYDVPDYA) from influenza hemagglutinin epitope coupled to KLH.
Specificity	This monoclonal antibody binds to HA-tagged recombinant protein.
Cross-reactivity	Not identified
Clonality	Monoclonal
Source	Mouse
Subclass	IgG3
Light Chain	Kappa
Formulation	Liquid, in pH7.4 PBS with 0.05% sodium azide.
Storage	Aliquot and store at -20°C for long term storage. Avoid repeated freeze and thaw cycles.
Applications	Western blotting (1:1,000). Immunoprecipitation and Immunostaining (1:200)
Research	Life science
Storage	Store at -20° C
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

Experimental data



Immunofluorescence staining of HA-tag fusion protein (transcription factor) in a stable expressing cell line (B) and control cell line (A).



HA Western Blot analysis of 293 cells transfected with HA-tagged vector (2) and untransfected control (1).

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited

S7.5(01)