



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Mouse anti-CA 50 Monoclonal Antibody Datasheet[Clone93]

Catalogue No	MO-T40014
Product name	Mouse anti human CA 50 monoclonal antibody [Clone 93]
Clone No	93
Product type	Primary antibodies.
Description	Mouse monoclonal antibody (MAb) to human CA 50, Protein G affinity purified
Target protein	Human CA 50.
Immunogen	Human CA 50 protein
Specificity	React with human CA 50.
Reactivity	Human, others not tested
Cross-reactivity	Not available
Clonality	Monoclonal
Source	Mouse
Subclass	IgG3; Kappa
Formulation	Lyophilized from a solution in 0.01M PBS, pH 7.2
Reconstitution	Double distilled water is recommended and to adjust the final concentration to 1.00 mg/ml; After reconstitution, aliquot and storage in freezer [-20 ⁰ C]. Avoid thaw and freeze cyclically
Applications	ELISA, antibody concentration 0.2 ng/ml can recognize coated 1 ug/ml of CA50/well.
Research	Oncology
Storage	Store at -20 ⁰ C
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

Lot
Expiration
Quantity

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited