



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Mouse anti-CA 125 Monoclonal Antibody Datasheet[Clone 23]

Catalogue No	MO-T40012C
Product name	Mouse anti CA 125 monoclonal antibody [Clone 23]
Clone No	23
Product type	Primary antibodies.
Description	Mouse monoclonal antibody (MAb) to human cancer antigen 125, Protein G affinity purified
Target protein	Human CA125.
Immunogen	Human CA 125.
Specificity	React with human CA 125.
Reactivity	Human, others not tested
Cross-reactivity	Not available
Clonality	Monoclonal
Source	Mouse
Subclass	IgG2a
Formulation	Lyophilized from a solution in 0.01M PBS, pH 7.2
Reconstitution	Double distilled water is recommended and to adjust the final concentration to 1.00 mg/ml
Applications	EIA
Research	Oncology
Storage	Store at -20° C
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited

S7.5(01)