



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Mouse anti-Human Free PSA Monoclonal Antibody-Tracer Datasheet

Catalogue No	MO-T40052A
Product name	Mouse anti-human PSA monoclonal antibody.
Clones	103
Product type	Matched antibody pair, tracer
Description	Mouse monoclonal antibody to human PSA, Protein G affinity purified
Target protein	Human Free PSA
Immunogen	Purified human PSA
Specificity	This monoclonal antibody reacts with both complexed and free human PSA
Reactivity	Human, others not tested
Cross-reactivity	Cross-react with complexed human PSA
Clonality	Monoclonal
Source	Mouse
Subclass	IgG2a
Light chain	Kappa
Formulation	Lyophilized from a solution in 0.01M PBS, pH 7.0
Reconstitution	Double distilled water is recommended to adjust the final concentration to 1.00 mg/ml
Applications	ELISA: This monoclonal antibody can be used as a tracer/detection antibody in a sandwich ELISA for human free PSA detection in combination with a free PSA specific capture antibody (Cat. No.: MO-T40052B).
Storage	Store at -20° C
Research	Oncology
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited

S7.5 (01)