



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Mouse anti-Human PMA Monoclonal Antibody [YPMA-2] Datasheet

Catalogue No	MO-T40083B
Product name	Mouse anti-human PMA monoclonal antibody [YPMA-2]
Clone No	YPMA-2
Product type	Primary antibodies
Description	Mouse monoclonal antibody (MAb) to human prostate mucin antigen (PMA), Protein G affinity purified
Target protein	Human prostate mucin antigen (PMA)
Immunogen	Membrane protein extracted from human prostate mucin tissues
Specificity	MAb YPMA-1 reacts specifically with prostate mucin antigen
Reactivity	Human, others not tested
Cross-reactivity	Does not react with fetal or benign prostate specimens, non-prostate carcinomas and normal human tissues
Clonality	Monoclonal
Source	Mouse
Subclass	IgG3
Formulation	Lyophilized from a solution in 0.01M PBS, pH 7.2
Reconstitution	Double distilled water is recommended and to adjust the final concentration to 1.00 mg/ml
Applications	Immunoperoxidase staining, immunofluorescent staining, and western blot
Research	Oncology and Cancer
Storage	Store at -20° C
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administration to humans and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited