



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Mouse anti-Human PSP Monoclonal Antibody [YPSP-1] Datasheet

Catalogue No	MO-T40085A
Product name	Mouse anti-human PSP monoclonal antibody [YPSP-1]
Clone No	YPSP-1
Product type	Primary antibodies
Description	Mouse monoclonal antibody (MAb) to human prostate secretory protein (PSP), Protein G affinity purified
Target protein	Human prostate secretory protein (PSP)
Immunogen	Highly purified PSP
Specificity	Reactive with PSP or human prostate epithelial cells
Reactivity	Human, others not tested
Cross-reactivity	Does not react with a variety of human normal tissues
Clonality	Monoclonal
Source	Mouse
Subclass	IgG1
Formulation	Lyophilized from a solution in 0.01M PBS, pH 7.2
Reconstitution	Double distilled water is recommended and to adjust the final concentration to 1.00 mg/ml
Applications	Immunoperoxidase staining, Western Blot and ELISA
Research	Oncology
Storage	Store at -20° C
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited

S7.5(01)