



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Mouse anti-CA 125 Monoclonal Antibody Datasheet[Clone110]

Catalogue No	MO-T40012A
Product name	Mouse anti CA 125 monoclonal antibody [Clone 110]
Clone No	110
Product type	Primary antibodies.
Description	Mouse monoclonal antibody (MAb) to human CA 125, Protein G affinity purified
Target protein	Human CA 125.
Immunogen	Human CA 125 protein
Specificity	React with human CA 125.
Reactivity	Human, others not tested
Cross-reactivity	Not available
Clonality	Monoclonal
Source	Mouse
Subclass	IgG1
Formulation	Lyophilized from a solution in 0.01M PBS, pH 7.2
Reconstitution	Double distilled water is recommended and to adjust the final concentration to 1.00 mg/ml
Applications	ELISA, antibody concentration 0.125 ng/ml can recognize coated 5 international unit of CA125/well.
Research	Oncology
Storage	Store at -20° C
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited

S7.5(01)