



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

HRP conjugated Goat anti-Human CRP Polyclonal Antibody Datasheet

Catalogue No	AI-20009P
Product name	HRP conjugated Rabbit anti-human C Reactive Protein (CRP) polyclonal antibody.
Product type	Primary antibody. As a tracer, matched with mouse anti-CRP monoclonal antibody [clone 63F4]
Description	Goat polyclonal antibody to human CRP, Protein G affinity purified, HRP conjugated.
Target protein	Human CRP
Immunogen	Human C Reactive Protein (CRP) from Pleural Fluid.
Specificity	This antibody recognizes human CRP.
Reactivity	Human, others not tested.
Clonality	Polyclonal
Source	Goat
Formulation	Liquid form. 0.01M PBS, pH 7.0 ± 0.1 with 50% glycerol
Applications	ELISA, Using as tracer with our matched mouse anti-CRP capturer antibody [Cat No: MO-M40023A],
Research	Biochemistry, inflammatory.
Storage	Store at -20° C. Aliquot and store at -20°C.
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited

Rev: S7.5(01)