



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

FITC Conjugated Goat anti-human IgG (Fc) Polyclonal Antibody Datasheet

Catalogue No	AI20001FI
Product name	FITC conjugated goat anti-human IgG (Fc) polyclonal antibody
Product type	FITC conjugated polyclonal antibody
Description	Produced by hyperimmunizing the host animal with highly purified human IgG (Fc) fraction. The goat IgG fraction is prepared from antisera with a Protein G affinity column. The polyclonal antibody is conjugated with FITC.
Target protein	Human IgG (Fc)
Immunogen	Highly purified human IgG (Fc) fraction
Specificity	This antibody recognizes human IgG (Fc) fraction
Reactivity	Human, others not tested
Clonality	Polyclonal
Source	Goat
Formulation	Liquid. Buffer: 10mM PBS, pH 7.2, 0.1% sodium azide as preservative
Applications	ELISA
Storage	Store at -20° C
Research	Immunology
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administration to humans and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited