



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Goat anti-Hepatitis B Surface Antigen (HBsAg) Polyclonal Antibody (Purified) Datasheet

Catalogue No	AI20002-P
Product name	Purified goat anti-human HBsAg Polyclonal antibody
Product type	Primary antibodies
Description	Goat anti-human hepatitis B Surface Antigen (HBsAg) polyclonal antibody, purified with a Protein G affinity column
Target protein	Human Hepatitis B surface antigen
Immunogen	Purified HBsAg from virulates
Specificity	This antibody recognizes Human Hepatitis B surface antigen
Reactivity	Human, others not tested
Clonality	Polyclonal
Source	Goat
Formulation	Lyophilized from a solution of 0.02M PBS, Ph7.4
Reconstitution	Use distilled water to reconstitute to desired concentration
Applications	ELISA
Storage	Store at -20° C
Research	Virology
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administration to humans and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited