



## **ANOGEN - A Division of YES Biotech Laboratories Ltd.**

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

### **Goat anti-HIV-1 Polyclonal Antibody (Purified) Datasheet**

<b>Catalogue No</b>	AI20003P
<b>Product name</b>	Purified goat anti-HIV-1 purified polyclonal antibody
<b>Product type</b>	Primary antibodies
<b>Description</b>	Goat anti-HIV-1 polyclonal antibody, purified with a Protein G affinity column
<b>Target protein</b>	HIV-1
<b>Immunogen</b>	Highly purified human immunodeficiency virus lysate
<b>Specificity</b>	This antibody recognizes: P17, P24, P45, P53, P55, P64, P120 and gp160. Strong positive bands on P17, P24 and P64
<b>Reactivity</b>	Human, others not tested
<b>Clonality</b>	Polyclonal
<b>Source</b>	Goat
<b>Formulation</b>	Lyophilized from a solution of 0.02M PBS, Ph7.4
<b>Reconstitution</b>	Use distilled water to reconstitute to desired concentration
<b>Applications</b>	ELISA & Western Blot
<b>Storage</b>	Store at -20° C
<b>Research</b>	Virology
<b>References</b>	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

**This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.**

**Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited**

S7.5(01)