



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

FITC Conjugated Goat anti-Mouse IgG Polyclonal Antibody Datasheet

Catalogue No	AI20004FI
Product name	FITC conjugated goat anti-mouse IgG polyclonal antibody
Product type	FITC conjugated polyclonal antibody
Description	Produced by hyperimmunizing the host animal with highly purified mouse IgG fraction. The goat IgG fraction is prepared from antisera with a Protein G affinity column. This polyclonal antibody is conjugated to FITC at a final FITC/protein ratio of 2:5.
Target protein	Mouse IgG
Immunogen	Highly purified mouse IgG fraction
Specificity	This antibody recognizes mouse IgG fraction
Reactivity	Mouse, others not tested
Clonality	Polyclonal
Source	Goat
Formulation	Liquid. Buffer: 10mM PBS, pH 7.2, 0.1% sodium azide as preservative
Applications	ELISA
Storage	Store at -20° C
Research	Immunology
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administration to humans and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited