



**ANOGEN** - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

**Goat anti-Mouse IgG Polyclonal Antibody (Purified) Datasheet**

<b>Catalogue No</b>	AI20004
<b>Product name</b>	Purified goat anti-mouse IgG polyclonal antibody
<b>Product type</b>	Polyclonal antibody
<b>Description</b>	Produced by hyperimmunizing the host animal with highly purified mouse IgG fraction. The goat IgG fraction is prepared with a Protein G affinity column
<b>Target protein</b>	Mouse IgG
<b>Immunogen</b>	Highly purified mouse IgG fraction
<b>Specificity</b>	This antibody recognizes mouse IgG fraction
<b>Reactivity</b>	Mouse, others not tested
<b>Clonality</b>	Polyclonal
<b>Source</b>	Goat
<b>Formulation</b>	Lyophilized from a solution of 0.02M PBS, pH7.4
<b>Reconstitution</b>	Use distilled water to dilute to desired concentration
<b>Applications</b>	ELISA
<b>Storage</b>	Store at -20° C
<b>Research</b>	Immunology
<b>References</b>	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

**This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administration to humans and animals for use in diagnostic and therapeutic procedures.**

**Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited**