



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Goat anti-Hepatitis C Virus (HCV) Polyclonal Antibody Datasheet

Catalogue No	AI20007
Product name	Goat anti-human HCV Polyclonal antibody
Product type	Primary antibodies
Description	Goat anti-human hepatitis C virus (HCV) polyclonal antibody is a antiserum obtained by immunizing the host with HCV polyprotein
Target protein	Human Hepatitis C virus
Immunogen	Recombinant Chimeric HCV polyprotein (60 Kd)
Specificity	This antibody recognizes: Recombinant Chimeric HCV Polyprotein (60Kd.) Recombinant CPC* + EPM** (1 a.a. -142 a.a.) Synthetic CPC* (1 a.a. -61 a.a.) Recombinant NS-3 (1252 a.a. -1477 a.a.) Synthetic NS-3 (1378 a.a.-1458 a.a) Synthetic NS-4a (1689a.a-1735a.a) * Capsid Protein C, ** Envelope protein M
Reactivity	Human, others not tested
Clonality	Polyclonal
Source	Goat
Formulation	Serum
Applications	ELISA & Western Blot (using 1:5000 dilution)
Storage	Store at -20° C
Research	Virology
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administration to humans and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited