



## **ANOGEN - A Division of YES Biotech Laboratories Ltd.**

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

### **Rabbit anti-Human Immunodeficiency Virus (HIV)-P24 Polyclonal Antibody Datasheet**

<b>Catalogue No</b>	AI20008
<b>Product name</b>	Rabbit anti-Human Immunodeficiency Virus (HIV)-P24 Polyclonal antibody
<b>Product type</b>	Primary antibodies
<b>Description</b>	Rabbit anti-HIV P24 polyclonal antibody is an anti-serum obtained by immunizing the host with recombinant HIV P24.
<b>Target protein</b>	HIV-P24
<b>Immunogen</b>	Highly purified recombinant human HIV-P24
<b>Specificity</b>	This antibody recognizes both native and recombinant HIV-P24.
<b>Reactivity</b>	Human, others not tested
<b>Clonality</b>	Polyclonal
<b>Source</b>	Rabbit
<b>Formulation</b>	Serum
<b>Applications</b>	ELISA
<b>Storage</b>	Store at -20° C
<b>Research</b>	Virology
<b>References</b>	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

**This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administration to humans and animals for use in diagnostic and therapeutic procedures.**

**Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited**

**S7.5(01)**