



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Rabbit anti-Human Immunodeficiency Virus (HIV)-P24 Polyclonal Antibody Datasheet

Catalogue No	AI20008P
Product name	Rabbit anti-Human Immunodeficiency Virus (HIV)-P24 Polyclonal antibody
Product type	Primary antibodies
Description	Rabbit anti-HIV P24 polyclonal antibody, purified with a Protein G affinity column
Target protein	HIV-P24
Immunogen	Highly purified recombinant human HIV-P24
Specificity	This antibody recognizes both native and recombinant HIV-P24.
Reactivity	Human, others not tested
Clonality	Polyclonal
Source	Rabbit
Formulation	Lyophilized from a solution of 0.02M PBS, pH 7.4
Applications	ELISA
Storage	Store at -20° C
Research	Virology
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administration to humans and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited

S7.5(01)