



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Rabbit anti-Mouse IgG Polyclonal Antibody Datasheet

| | |
|-----------------------|--|
| Catalogue No | AI50005 |
| Product name | Rabbit anti-mouse IgG Polyclonal antibody |
| Product type | Primary antibodies |
| Description | Produced by hyperimmunizing the host animal with highly purified mouse IgG fraction. The rabbit IgG fraction is prepared from mono-specific antisera (NH ₄) ₂ SO ₄ precipitation followed by ion exchange chromatography |
| Target protein | mouse IgG |
| Immunogen | Highly purified mouse IgG fraction |
| Specificity | This antibody recognizes mouse IgG |
| Reactivity | Mouse, others not tested |
| Clonality | Polyclonal |
| Source | Rabbit |
| Formulation | Lyophilized from a solution of 0.02M PBS, pH7.4 |
| Applications | ELISA |
| Storage | Store at -20° C |
| Research | Immunology |
| References | If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude. |

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administration to humans and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited