



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Rabbit anti-Human Growth Hormone (hGH) Polyclonal Antibody Datasheet

Catalogue No	AI50006P
Product name	Purified rabbit anti-human growth hormone polyclonal antibody
Product type	Primary antibodies
Description	Rabbit anti-human growth hormone polyclonal antibody, purified with a Protein G affinity column
Target protein	hGH
Immunogen	Recombinant Human growth hormone purified by HPLC. The purity of rhGH was 98% as determined by SDS-PAGE
Specificity	This antibody recognizes human growth hormone
Reactivity	Human, others not tested
Clonality	Polyclonal
Source	Rabbit
Formulation	Lyophilized from a solution of 0.01M PBS, pH 7.4
Reconstitution	Use double distilled water to reconstitute to desired concentration
Applications	ELISA
Storage	Store at -20° C
Research	Immunology
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administration to humans and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited