



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Rabbit Anti Phospho Bcl-2 (Thr56) Polyclonal Antibody Datasheet

Catalogue No	AI70003B
Product name	Purified rabbit anti phospho Bcl-2 (Thr56) polyclonal antibody
Product type	Primary antibody
Description	Rabbit anti phospho Bcl-2 (Thr56) polyclonal antibody
Target protein	Thr56 phosphorylated B cell lymphoma-2 (Bcl-2).
Immunogen	A synthetic peptide (SQPGHT _p PHPASR) derived from Bcl-2 phosphorylated at threonine 56, conjugated with KLH for immunization.
Specificity	This antibody recognizes the Thr56 phosphorylated Bcl-2 peptide.
Reactivity	Thr56 phosphorylated Bcl-2 peptide of human and similar molecule of other species.
Clonality	Polyclonal
Source	Rabbit
Formulation	Lyophilized in pH 7.4, 0.01M PBS. To reconstitute, add appropriate volume of dH ₂ O to reach a final antibody concentration of 1mg/ml.
Applications	ELISA: dilute 1:2,000 – 1:100,000 to detect the Thr56 phosphorylated Bcl-2 peptide coated at 100ng/well. Other applications have not evaluated.
Storage	Store in aliquots at -20°C. Avoid repeated freeze/thaw cycles
Research	Apoptosis signaling pathway
References	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited