




**ANOGEN** - A Division of YES Biotech Laboratories Ltd.  
2355 Derry Road East, Unit 23, Mississauga, ON L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

### Rabbit anti-Human IgA Polyclonal Antibody Datasheet

<b>Catalogue No</b>	AI80004
<b>Product name</b>	Purified rabbit anti-human IgA polyclonal antibody
<b>Product type</b>	Primary antibodies
<b>Description</b>	Rabbit anti-human IgA polyclonal antibody, affinity purified.
<b>Target protein</b>	Human IgA
<b>Immunogen</b>	Human IgA
<b>Specificity</b>	This antibody mainly recognizes human IgA heavy chain.
<b>Reactivity</b>	Human, others not tested
<b>Clonality</b>	Polyclonal
<b>Source</b>	Rabbit
<b>Formulation</b>	Liquid, in 0.01M PBS, pH7.2 containing 0.01% BSA, 0.05% sodium azide and 50% glycerol
<b>Applications</b>	Suggested application and working condition: Western blot: 1:80,000
<b>Storage</b>	Store in aliquots at -20°C. Avoid repeated freeze/thaw cycles
<b>Research</b>	Biochemistry.
<b>References</b>	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide any product of choice as gratitude.

### Western blot:



**This product is for LABORATORY RESEARCH USE and further manufacture only, and can not be administrated to human and animals for use in diagnostic and therapeutic procedures.**

**Manufactured by ANOGEN - A Division of YES Biotech Laboratories Limited**